



2-CAN™



UNIQUE BIOLOGICS DELIVERY

to Virtually any Injured Bone

A disposable answer to the Biologics Era, 2-Can, our Multi-Channeled Side-Port Cannula, fulfills your need for biologics delivery.

- Passing over a guide wire, the **multi-channeled design isolates the biologics delivery chamber** from the guide wire preventing unwanted biologic extrusion while keeping the guide wire as clean as possible.
- Creating new treatment paradigms for problem fractures, **2-Can is available in a range of sizes** ensuring the success of almost every procedure.
- An integral plunger **maximizes the efficiency of biologic delivery** and assists in maintaining an open implant channel.
- Use 2-Can after drilling/tapping, but prior to screw or nail placement. Precision delivery is achieved with **2-Can's Side-Port Cannula design** with a radiopaque tip.

Flow-FX, LLC - 815.531.4424 - 9301 West 191st Street, Mokena, IL 60448 - www.Flow-FX.net





2-Can™ 8L
103-1-04 - 8mm O.D. x 450mm



2-Can™ 8S
103-1-05 - 8mm O.D. x 350mm



2-Can™ 6.5
103-1-03 - 4.7mm O.D. x 220mm



2-Can™ 5.5
103-1-02 - 3.5mm O.D. x 170mm

