



FLOW-NAIL™ WITH LAGLOCK™ INFORMATION SHEET



INTRAMEDULLARY FIXATION REVOLUTIONIZED WITH PRECISION DELIVERY OF BONE VOID FILLER AND SLIP-PROOF, STATIC LOCKING LAGLOCK FEATURE



Patented

Patented, 510k-cleared Flow-Nail offers exclusive advantages:

- Patented side port cannula delivers bone void filler through fenestrations while protecting the hip joint from unintended intrusion of biologic material
- Advanced delivery mechanism ensures precise placement of bone void filler to the implantation site using manual low-pressure delivery
- Invented for fracture fixation, Flow-Nail gives the surgeon an option regarding standalone construct or supplementation with bone void filler
- Available in 125° and 130° neck angles with 180mm and 220mm lengths to allow for the treatment of the majority of fractures as a standalone device
- The Laglock lag screw creates a unique, patent-pending, static nail construct for challenging fractures. It withstands a 400-pound static load in line with the screw and 1 million cycles at 72nm fatigue tests.



Flow-Nail precisely delivers bone void filler where you need it while protecting the joint from unintended intrusion of biologic material.

www.Flow-FX.net



REIMAGINING THE SCIENCE OF SKELETAL HEALTH



Trochanteric Nail

Catalog #	Description
101NAIL-125-18011	11mm/125°/180mm
101NAIL-125-22011	11mm/125°/220mm
101NAIL-130-18011	11mm/130°/180mm
101NAIL-130-22011	11mm/130°/220mm
101NAIL-125-18010	10mm/125°/180mm
101NAIL-130-18010	10mm/130°/180mm
101NAIL-125-22010	10mm/125°/220mm
101NAIL-130-22010	10mm/130°/220mm



Cortical Screw

Catalog #	Description
C20-C80	20-60mm (2mm Increments) 60-80mm (5mm Increments)

Plunger

Catalog #	Description
102-51-02D	6.5/7.3 Single use Plunger



PEG

Catalog #	Description
PEG	Collapse Limiter Peg

Cap

Catalog #
101CAP-16-01TI



Lag Screw with Laglock

Catalog #	Description
101LGLK-1080	80mm Lag Lock Screw
101LGLK-1085	85mm Lag Lock Screw
101LGLK-1090	90mm Lag Lock Screw
101LGLK-1095	95mm Lag Lock Screw
101LGLK-1100	100mm Lag Lock Screw
101LGLK-1105	105mm Lag Lock Screw
101LGLK-1110	110mm Lag Lock Screw
101LGLK-1115	115mm Lag Lock Screw
101LGLK-1120	120mm Lag Lock Screw



Lag Screw (non-locking)

Catalog #	Description
101LSCR-(1080TI-1120TI)	80-120mm (5mm Increments)



Anti-Rotation Screw

Catalog #	Description
A50-100	50-100mm (10mm Increments)



Find out how Flow-FX products can help improve patient outcomes and lower costs in your practice. Visit Flow-FX.net to learn more.

www.Flow-FX.net